

Product Name : Surgery Microscope	Product Code : NLABMEDICAL-M3300002
Description : Surgery Microscopes India, Surgery Microscope manufacturer India, suppliers, Medical Surgery Microscopes India, Hospital Surgery Microscope India	
Technical Specification : Reliability Dimensionally accurate Accurate testing results Innovative designed, highly versatile with latest 45° inclined optical system 15x wide field eye pieces 3 Step magnification (5x, 10x & 20x) or 5 Step magnification (0. 4x, 0. 63x, 1x, 1. 6x & 2. 5x) or Continuous zoom (5x ~ 25x) Beam splitter with 'c' mount CCD digital camera Working distance: f= 200mm objective lens Interpupillary distance adjustment: 55mm to 75mm Focusing – manual or motorized with foot control LCD monitor/screen: mount on microscope to display live surgery Light source : 15v-150w / 24v-250w twin halogen reflector lamp	

Illumination control : continuous adjustment

Co-axial through-the-lens fiber optic light guide, white bright light

Field of illumination : 70mm

Filter : cobalt blue, red-free, heat absorbing filter

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003