

**Product Name :**  
Surgical Cautery Unit

**Product Code :**  
NLABMEDICAL-M3300001



**Description :**

Surgical Cautery Unit

**Technical Specification :**

Monopolar cutting, coagulation, blended output

Bipolar cutting / Bipolar coagulation

Independent activation of hand switch of foot switch

Audio visual alarm for separate mode

Special PVC foot pedal for higher corrosion resistance

Monopolar Generator:

Frequency: 400+ - 25 Khz

Cut Control (Pure): 250 Watts 500 Ohm Load

Cut Control (Blend): 175 Watts 500 Ohm Load

Coagulation Control: 160 Watts 500 Ohm Load

Bipolar Generator:

Frequency: 650 + - 25 Khz

Bipolar Coagulation: 80 Watts 125 Ohm Load

Dimensions: W 285 mm X D 295 mm X H 120 mm

Power Consumption: 750 VA

Gross Weight: 20 lbs (approx)

Voltage Required: 220V AV 50 Hz

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003