

<b>Product Name :</b> Syringe Infusion Pump	<b>Product Code :</b> NLABMEDICAL-M3320002
<b>Description :</b>	
Syringe Infusion Pump	
<b>Technical Specification :</b>	
Self-explanatory with abbreviated operating instructions on LCD.	
Easy-to-read liquid-crystal displays all essential data in large size numbers.	
Reducing the chance of erroneous readings.	
LCD provided with back light for viewing in dark environments.	
Selectable three-step occlusion detection pressures.	
Low, mid & high-occlusion setting.	
"Keep-Vein-Open" function.	
Infusion rate in KVO is 1% but not less than 0. 1 ml/hr.	
This prevent the infusion line from blocking due to blood clots.	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003