

Product Name :
Syringe Pump Dual Channel

Product Code :
NLABMEDICAL-M3420001



Description :

Syringe Pump Dual Channel

Technical Specification :

Syringe Pump Dual Channel

Dual-channel with different flow rate

Channel sleeping function

3-level occlusion alarm limits setting

Limits setting

Automatically reduce impact dose release after occlusion

Delivered volume inquiry

KVO flow rate

Fast rate control

Purge button activate safe protection

Flow rate setting button lock

Abundant optic and sound alarm

AC and internal power indication

Standard RS232 interface

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003