

Product Name :
Table Top Front Loading Autoclave

Product Code :
NLAB-LABORATORYEM47010



Description :

Table Top Front Loading Autoclave

Technical Specification :

Table Top Front Loading Autoclave -

Features:-

Table Top Front loading Sterilizer / Autoclave.

Digital Pressure Display & also Analogue pressure display.

Auto drain of water to reservoir tank after cycle completion.

Completely automatic & with error message.

Seamless chamber.

Feasibility to take out instruments from chamber even after electricity fails.

Fast cycle time with Digital timer for Wet & Dry Cycle.

Mode Select:-

Instrument programs to sterilize metallic devices sterilization temperature : 134°C, sterilization time : 4 minutes & dry time : 15 minutes.

Liquid programs to sterilize liquid sterilization temperature : 121°C, sterilization time : 20 minutes & dry time : 15 minutes.

Rubber programs to sterilize rubber & plastic sterilization temperature : 121°C, sterilization time : 20 minutes & dry time : 15 minutes.

Manual Program:-

Manual Program Manual program as per user requirements program can set pre vacuum, post vacuum sterilization temperature : 121°C~134°C, sterilization time : 1~60 minutes, dry time : 1~60 minutes.

We are leading manufacturers, suppliers of Table Top Front Loading Autoclave for Autoclave. Contact us to get high quality Table Top Front Loading Autoclave for Autoclave for schools, colleges, universities, research labs, laboratories and various industries.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003