

Product Name :
Tele Thermometer

Product Code :
NLAB-SCIENTIFICEQ28026



Description :

Tele Thermometer

Technical Specification :

Use for Pyrogen testing 0 to 6 channel (Digital display Temperature) with 6 probes.

Range: 0 to 60°C

Resolution: 0.1°C

Accuracy : $\pm 0.2\%$ F.S ± 1 Digit

Input: RTD (PT-100) or IC Sensor

Channels: 6 Channels Switch Selectable

Display: 3 1/2 digit red LED display

Power: 220V $\pm 10\%$ AC 50 HZ

Dimension: 275x175x75 mm

Weight: 2 Kgs

Accessories:

6 Pyrogen Probes

Dust Cover

Instruction Manual

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003