

**Product Name :**  
 Telephone Trainer

**Product Code :**  
 COM0001



**Description :**

Telephone Trainer

**Technical Specification :**

Telephone Trainer For Electrical Lab Training.

Telephone Trainer

The

Telephone Trainer presenting by Vinytics is a complete kit to understand the working of telephone instruments with the idea of working of exchange.

The trainer has DIP switches to create the faults in the instruments and test points to give the complete idea about the servicing/working of telephone instrument.

Based on Basic Telephone Instrument Circuits.

16 Fault creating switches for different sections.

Circuit

Diagram on board describing different sections in telephone like Line protector, Voltage dropper, Ringer, Dialer, Key matrix, Sound amplifier section etc.

Test points to measure the signals in different sections.

Can be use as normal telephone instrument.

User's manual with fault finding

## **Naugralabequipments**

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-0171-2643080, 2601773