

Product Name : Temperature Controlled X -Ray Film Processing Unit	Product Code : NLABMEDICAL-M3350006
Description :	
Temperature Controlled X -Ray Film Processing Unit	
Technical Specification :	
Temperature Controlled X -Ray Film Processing Unit	
Available in request dimensions for capacity.	
Temperature maintained at 20 degree Celsius.	
The unit is made up of Stainless steel 20 swg.	
Od 316 grade.	
The unit is incorporated with suitable capacity sealed compressor and necessary accessories for maintaining 20 degree Celsius.	
The unit is having two compartments, one to maintain solution temperature at 20 degree Celsius and the second for developer, rinser, fixer and final wash for x-ray films with a provision of inlet, outlet & running water system.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003