

Product Name : Temperature Measurement and Calibration	Product Code : NLAB-TECHNICALAB36006
Description :	
Temperature Measurement and Calibration	
Technical Specification :	
<p>A bench top unit comprising a hypsometer /hot water bath and ice flask which generate accurate fixed points and variable temperatures</p> <p>An accurate platinum resistance thermometer with five point calibration certificate and temperature indicated directly in °C is included for reference.</p> <p>Key Features:</p> <p>Thermometric properties (physical changes due to temperature) of different sensors are investigated</p> <p>Condensing vapour (using a hypsometer) provides an accurate fixed point for calibration</p> <p>Stirred hot water bath provides variable temperatures</p> <p>Reference with calibration supplied</p> <p>All electrical sensors can be logged using a PC (optional teaching software available)</p>	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003