

Product Name :
Temperature Measurement and Calibration

Product Code :
NLAB-TECHNICALAB36006

Description :

Temperature Measurement and Calibration

Technical Specification :

A bench top unit comprising a hypsometer /hot water bath and ice flask which generate accurate fixed points and variable temperatures

An accurate platinum resistance thermometer with five point calibration certificate and temperature indicated directly in °C is included for reference.

Key Features:

Thermometric properties (physical changes due to temperature) of different sensors are investigated

Condensing vapour (using a hypsometer) provides an accurate fixed point for calibration

Stirred hot water bath provides variable temperatures

Reference with calibration supplied

All electrical sensors can be logged using a PC (optional teaching software available)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003