


Product Name :
Temperature Sensor

Product Code :
NLAB-VOCATIONALTVET28010



Description :

Temperature Sensor

Technical Specification :

$\pm 2^{\circ}\text{C}$ accuracy over temperature

Voltage Input: 2.7 V to 5.5 VDC

10 mV/ $^{\circ}\text{C}$ scale factor

Operating Range: -40°C to $+125^{\circ}\text{C}$

$\pm 0.5^{\circ}\text{C}$ linearity

We are well-known manufacturers, OEM suppliers of Temperature Sensor for Data Acquisition System. Contact us for high quality Temperature Sensor for Data Acquisition System for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com