

Product Name :
Tenoning Machine

Product Code :
NLAB-MACHINELABWK33029



Description :

Tenoning Machine

Technical Specification :

Power: 3 x 2.2kW (3 phase)

Main cutter heads: 2175 x 60mm; 1125 x 60mm

Vertical adjustment bottom block: 0-60mm

Vertical adjustment scribing head: 95mm

Vertical adjustment top block: 0 – 100mm

Horizontal adjustment scribing head: 95mm

Horizontal adjustment top block: +/- 32mm

Scribing head shaft: 75 x 30mm

Table size : 800 x 310mm

Come complete with knives (60mm x 12mm x 1.5mm – 10

And scribing head (125 x 65mm)

Wood Workshop Portable Tools Tenoning Machine Manufacturer, Supplier and Tenoning Machine Exporter in India

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003