

**Product Name :**  
Tensile Strength Testing Machine

**Product Code :**  
NLAB-MECHANICALAB250008



**Description :**

Tensile Strength Testing Machine

**Technical Specification :**

Tensile Strength Tester and elongation are the two prime characteristics of most of the raw materials, whether they are metals or non-metals such as foam, rubber, textiles, plastic, paper or finished products such as rods, wires, ropes, yarns, belts etcetera. These properties very often play the major role in determining the suitability of any raw material for any specific application. It is, therefore, of utmost importance to determine these characteristics accurately, conveniently, and quickly.

**Specifications:**

Range: 250 Kg with 0.1 kg least count digital display

Elongations: Up to 900 mm X 0.1 mm least count digital display

Motor: ½ Hp, 1425 rpm, 220 V AC single phase induction motor

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003