


Product Name :
Tensile Testing Equipments

Product Code :
NLAB-ENGINEERINGLB31005



Description :

Tensile Testing Equipments

Technical Specification :

Tensile Testing Machine Lab Equipments Manufacturer, Suppliers & Exporters Bench Top Tensile Testing Machine

Experiments

Tensile tests up to 20 kN on specimens made of different metals, to find material characteristics such as upper and lower yield strengths, tensile strength and overall extension.

Tests of Young's Modulus for the specimen material .

Operating Conditions

Operating environment: Laboratory

Storage temperature range: °C to +55°C (when packed for transport)

Operating temperature range: +5°C to +40°C

Operating relative humidity range: 80% at temperatures < 31°C decreasing linearly to 50% at 40°C

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003