

**Product Name :**  
Tensile Testing Equipments

**Product Code :**  
NLAB-ENGINEERINGLB31005



**Description :**

Tensile Testing Equipments

**Technical Specification :**

Tensile Testing Machine Lab Equipments Manufacturer, Suppliers & Exporters Bench Top Tensile Testing Machine

**Experiments**

Tensile tests up to 20 kN on specimens made of different metals, to find material characteristics such as upper and lower yield strengths, tensile strength and overall extension.

Tests of Young's Modulus for the specimen material .

**Operating Conditions**

Operating environment: Laboratory

Storage temperature range: °C to +55°C (when packed for transport)

Operating temperature range: +5°C to +40°C

Operating relative humidity range: 80% at temperatures < 31°C decreasing linearly to 50% at 40°C

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003