

**Product Name :**  
Tensile Testing Machine

**Product Code :**  
NLAB-ENGINEERINGLB31016



**Description :**

Tensile Testing Machine

**Technical Specification :**

High accuracy of  $\pm 1\%$

Three measuring ranges for accurate testing.

Easy in operation and maintenance.

Incorporate high level technology.

Highly Accurate and precise.

Modular design

Affordable

Safe and reliable

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003