

Product Name :
Ternary Systems Study Bench

Product Code :
NLAB-ENGINEERINGLB36025



Description :

Ternary Systems Study Bench

Technical Specification :

Thermostatically controlled water bath with the circulation pump.

3 cells manufactured from borosilicate glass with jacket and stirrer.

We are well-known manufacturers, OEM suppliers of Ternary Systems Study Bench for Chemical Engineering Lab Equipments. Contact us for high quality Ternary Systems Study Bench for Chemical Engineering Lab Equipments for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003