

Product Name :
Test Bench for Different Pumps

Product Code :
NLAB-EDUCATIONTVET14013



Description :

Test Bench for Different Pumps

Technical Specification :

Autonomous mobile unit for determining the characteristics of 8 different pumps (depending on accessories) at different speeds

5 pump locations (4 can be active at the same time)

Integrated control valve to study the effect of a loss on suction

API control panel for display of suction, flow, water height, speed and torque

Electronic measurement of flow and suction and discharge pressures (except retroactive pumps)

Supplied as standard with centrifugal pump (max. 137l / min, water height 9m) and gear pump (max. 13.7 l / min, water height 60m)

Provided as standard for data acquisition by PC

Recent configuration computer, supplied with 500VA inverter

Accessories:

Additional centrifugal pump, type C3-MKII-20SP (or similar),

Axial pump, type C3-MKII-22, max. 100 l / min, water height 1.8m

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003