

Product Name :
Test Duriez Normal Ø80 mm

Product Code :
NLAB-EDUCATIONTVET11011



Description :

Test Duriez Normal Ø80 mm

Technical Specification :

Mold body

Demoulding base

Upper and lower piston

Demoulding piston

Cleat set

For cold tests:

Grooved upper and lower piston

Grooved release piston

For the Duriez dilated Ø120 mm test:

Mold body

Demoulding base

Upper and lower piston

Demoulding piston

Cleat set

Grooved upper and lower piston

Grooved release piston

For the conservation of Duriez test specimens:

Climatic chamber

Stainless steel construction, interior glass door

Peltier technology

Temperature: 0 to 70 ° C

Humidity: 10 to 90% RH

Capacity: 256 liters

PID Twin Display controller, with 2 TFT color screens

Grid / tray runners: 9

Max load : 20kg per grid / tray - 200kg total

Internal dimensions (WxHxD): 64 x 80 x 50 cm

External dimensions (WxHxD): 82.4 x 118.3 x 75.6 cm

Power supply: 230V - 50Hz - 920W

Weight: 122 kg

Supplied with 2 stainless steel grids, factory calibration certificate at 10 ° C, 37 ° C and 30 ° C / 60% RH and water tank

Additional grid

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

