

Product Name :
Test Tube Borosilicate Glass

Product Code :
NLAB-EDUSCHOOLEQ46012



Description :

Test Tube Borosilicate Glass

Technical Specification :

We have a wide range of glass test tube in various capacities.

The test tubes are generally used in chemistry labs to hold, mix and heat small quantities of aqueous solutions.

These are made up of high quality glass to withstand different chemical reactions taking place in different experiments.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003