

Product Name :
Test Tube Brush W. Wool

Product Code :
NLAB-EDUSCHOOLEQ207005



Description :

Test Tube Brush W. Wool

Technical Specification :

Test tube brush with wool tip.

Twisted brush

Head of brush made of wool

Diameter: 20 mm

Handle: twisted eye

Total length: 270 mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

