

Product Name :
Test Tube Holder Wire Type

Product Code :
NLAB-EDUSCHOOLEQ46013



Description :

Test Tube Holder (Wire Type)

Technical Specification :

Catering to the requirements of our clients, we also provide wire type test tube holders.

These are made up of plated steel wire holder with finger grip and wire collar.

The alignment is done with precision so as to maintain a firm grip to the test tube.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003