

**Product Name :**  
Test Tube Holder Wooden

**Product Code :**  
NLAB-EDUSCHOOLEQ46014



**Description :**

Test Tube Holder (Wooden)

**Technical Specification :**

These test tubes are made up of wood and are well finished too.

Peg type test tube holders come with spring closing action using metal wire.

Our test tube holders are known for their durability, accurate designs and many other features.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003