

**Product Name :**  
Thermistor Characteristics Trainer

**Product Code :**  
NLAB-ELECTRICALAB210045



**Description :**

Thermistor Characteristics Trainer

**Technical Specification :**

Digital Meters:

Voltmeter 20VDC

Ammeter 200mA DC

Power supplies:

DC Supply IC Regulated 0-6V DC, 150mA.

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Fuse for Short Circuit protection

Components are mounted on the panels are:

Thermistor 10K.

Heating Oven On Board Controlled By ON/OFF Switch.

Voltage Control through Potentiometer.

110 Degree Thermometer

Salient Features:

Instruction manual.

Operated on Mains power 230V, 50Hz.

Patch Cords 4mm.

Supply provided through the 3 pin computer socket with appropriate fuse

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x8"

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003