

Product Name :
Thermistor Characteristics Trainer

Product Code :
NLAB-ELECTRICALLAB210045



Description :

Thermistor Characteristics Trainer

Technical Specification :

Digital Meters:

Voltmeter 20VDC

Ammeter 200mA DC

Power supplies:

DC Supply IC Regulated 0-6V DC, 150mA.

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Fuse for Short Circuit protection

Components are mounted on the panels are:

Thermistor 10K.

Heating Oven On Board Controlled By ON/OFF Switch.

Voltage Control through Potentiometer.

110 Degree Thermometer

Salient Features:

Instruction manual.

Operated on Mains power 230V, 50Hz.

Patch Cords 4mm.

Supply provided through the 3 pin computer socket with appropriate fuse

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x8"

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003