

Product Name :
Thermo Forming - Membrane Press

Product Code :
NLAB-MACHINELABWK15002



Description :

Thermo Forming - Membrane Press

Technical Specification :

Can be used for PVC, PET, PE & Mica

Overall size of Tray

Maximum Work Area for Job

Maximum thickness of Job

Maximum Depth of Tray

Pneumatic pressure required

Vasuum Pump mounted on machine

Vacuum Motor power

No. of Cylinders

Cylinder Stroke

Oil Pump motor power

Type of heating

Heating oil temperature

Heater Power

Overheat thermal cutout

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003