

Product Name :
Thermocouple Calibration Test Rig

Product Code :
NLAB-ENGINEERINGLB52012



Description :

Thermocouple Calibration Test Rig

Technical Specification :

DESCRIPTION

Here simple demonstration setup is provided to experiment with characteristics of Thermocouple. Here using this setup one can calibrate different types of thermocouples.

SPECIFICATIONS

1. Water Bath with heating arrangement.
2. Thermocouple with necessary signal conditioning circuit.
3. Thermometer.
4. Measurement of oC and Mili Volt on Digital Display
5. Digital Indicator: -

- Least Count: - 1 deg

- Digit: - 3 1/2

6. Powder coated enclosure.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003