

Product Name :
Thermostatic Water Baths (Stretching Inner Chamber)

Product Code :
NLAB-SCIENTIFICEQ34004



Description :

Thermostatic Water Baths (Stretching Inner Chamber)

Technical Specification :

Thermostatic Water Baths (Stretching Inner Chamber)

Used for precision constant temperature and auxiliary.

Characteristics:

Stainless steel stretching inner container and roof cover.

Intelligent temperature controller with PID program, digital tube display and function of timing and over-temperature protection.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003