

Product Name :
Thermostatic Water Baths

Product Code :
NLAB-SCIENTIFICEQ34005



Description :

Thermostatic Water Baths

Technical Specification :

Thermostatic Water Baths

Used for precision constant temperature and auxiliary

Characteristics:

two holes water bath, can be used as a digital two holes water bath.

six holes water bath, can be used as a digital six holes water bath.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003