

Product Name : Three Way Syringe	Product Code : NLABMEDICAL-M32100012
Description :	
Three Way Syringe India, Three Way Syringe manufacturer India, suppliers, Medical Three Way Syringe India, Hospital Three Way Syringe India	
Technical Specification :	
Offered come in quality construction finish and come with support of directly interchangeable finish with metal tips.	
Further, these come constructed with outer clear tube as well as with color inner.	
Some of its features include no Syringe conversion required before use, plastic material finish, more convenient for sterilization and cleaning, equipped with a special anti-siphon device, tube is not susceptible to plugging.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003