

**Product Name :**  
Tilting Flume

**Product Code :**  
DM541



**Description :**

Tilting Flume

**Technical Specification :**

A 76mm wide, 250mm high open channel.  
Clear acrylic sides to give visibility of the working section.  
A jacking system permits the slope of the channel bed to be adjusted between -1% and +3%.  
Service module: Made from Acrylic  
Sump tank capacity: 250 litres (S. S.)  
Volumetric tank capacity: 40 litres (GRP)  
Inlet tank with flow stilling arrangement. (GRP)  
Include a Venturi flume, sharp and broad crested weirs, 2 vernier level gauges, adjustable undershot weir and crump weir.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003