

Product Name :
Tilting Flume

Product Code :
DM541



Description :

Tilting Flume

Technical Specification :

A 76mm wide, 250mm high open channel.
Clear acrylic sides to give visibility of the working section.
A jacking system permits the slope of the channel bed to be adjusted between -1% and +3%.
Service module: Made from Acrylic
Sump tank capacity: 250 litres (S. S.)
Volumetric tank capacity: 40 litres (GRP)
Inlet tank with flow stilling arrangement. (GRP)
Include a Venturi flume, sharp and broad crested weirs, 2 vernier level gauges, adjustable undershot weir and crump weir.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003