

**Product Name :**  
Time Concepts Instruments For Maths Lab

**Product Code :**  
NLAB-EDUSCHOOLEQ169040



**Description :**

Time Concepts Instruments For Maths Lab

**Technical Specification :**

**Dummy Clock**

This Fiber dummy clock is used to teach time concepts with moveable hand of hours & minuets. Made up to light weight fiber sheet, easy to handle.

**Student Clock Write And Wipe**

Encourage class to participation with write-on/wipe-off clock. This set of 5 clocks is Great for small-group activities. clocks with movable hands and a place to write in digital time. Great way to encourage student participation in time-writing activities.

**Geared Teacher Clock**

This Geared clock Accurate to the hours with the movement of minute Arm. This geared demonstration clock is a great manipulative of time-telling skills for your students. This clock is provided with bright color clock arms, featuring easy-to-read. Hidden gears maintain correct hour and minute relationships as you manipulate movable

hands. Clock is made of durable plastic.

#### Student Time Indicator

To compare local times of different countries with GMT. Read different times for different countries and compare.

#### Time & Work Kit

This kit is provided with assorted color right angle triangle with digital stop watch. This kit help to calculate work dine in a same time span by the different group.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003