

Product Name :
Time Concepts Instruments For Maths Lab

Product Code :
NLAB-EDUSCHOOLEQ169040



Description :

Time Concepts Instruments For Maths Lab

Technical Specification :

Dummy Clock

This Fiber dummy clock is used to teach time concepts with moveable hand of hours & minuets. Made up to light weight fiber sheet, easy to handle.

Student Clock Write And Wipe

Encourage class to participation with write-on/wipe-off clock. This set of 5 clocks is Great for small-group activities. clocks with movable hands and a place to write in digital time. Great way to encourage student participation in time-writing activities.

Geared Teacher Clock

This Geared clock Accurate to the hours with the movement of minute Arm. This geared demonstration clock is a great manipulative of time-telling skills for your students. This clock is provided with bright color clock arms, featuring easy-to-read. Hidden gears maintain correct hour and minute relationships as you manipulate movable

hands. Clock is made of durable plastic.

Student Time Indicator

To compare local times of different countries with GMT. Read different times for different countries and compare.

Time & Work Kit

This kit is provided with assorted color right angle triangle with digital stop watch. This kit help to calculate work dine in a same time span by the different group.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003