

Product Name :
Tincture Press

Product Code :
NLAB-SCIENTIFICEQ28027



Description :

Tincture Press

Technical Specification :

Extra heavy duty for making tinctures, decoctious infusion etc. cap. 1 Lit. stainless steel disc, with screw drive, stainless steel mesh as inner perforated vessel with stop cock.

We further offer a range of Tincture Press, which is conceptualized for separation of marc of herbs and extracts made by percolation / maceration / extraction of herbs & solvents by pressing method. The unit consists of double chamber that includes a plunger pressing plate of approx. 1/2 inch of thickness, supposed to move inside the chamber with the support of hand operated mechanical rotation pressing system. The unit is further attached with an outlet valve. The whole assembly is mounted on heaving M.S. base stand with M.S. / S.S. sieve for filtration purpose.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003