

Product Name :
Tissue Embedding System

Product Code :
NLAB-LABORATORYEM28034



Description :

Tissue Embedding System

Technical Specification :

Setting temperature range : ambient to 99 deg C.

Five independent heating zones for better performance.

Programmable for both ON TIME & OFF TIME.

Electrically heated forceps holder.

Set temperature memorization.

Panel temp. Display for all heated / cooled areas.

Paraffin flow control can be done manually by a switch or through foot switch.

Removable waster tray.

Bright illumination of the work area.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003