

Product Name :
Tissue Flotation Bath

Product Code :
NLAB-EDUSCHOOLEQ29013



Description :

Tissue Flotation Bath

Technical Specification :

Tissue Flotation Bath is useful to float tissue wax sections to flatten for staining process.

Double Walled chamber, inner made of black anodised aluminum and outer wall made of mild steel sheet duly painted.

Gap between two insulated with glass wool.

Temperature ambient to 700° C 20° C controlled by a thermostate.

It is supplied completely with Pilot Lamp, Temperature Control Knob, Cord-Plug etc.

To work on 220V AC 50Hz single phase.

Chamber Size:- Diameter 225 mm, Depth 70 mm, Rim 40 mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003