

**Product Name :**  
Toe Cap Impact Tester Machines

**Product Code :**  
NLAB-TESTINGLABM57001



**Description :**

Toe Cap Impact Tester Machines

**Technical Specification :**

Toe Cap Impact Tester is made on a horizontal platform which has two guide bars which guide the impact weight. There is a mild steel plunger of specified dimensions which is used the intermediate parts between the impact weight and the toe cap. The plunger is freely supported by the toe cap. The base of the machine is placed on the hard wood of specified dimensions. A capsule for measuring the clearance at the time of maximum depression is also supplied with the instrument. All the parts are either painted or chrome plated to give the instrument a corrosion free long life

**Specifications:**

Clamping Block: 150 x 150 x 19 mm, hardness at least 60 HRC

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

