

Product Name :
Toggle Press

Product Code :
NLAB-MACHINELABWK26008



Description :

Toggle Press

Technical Specification :

Toggle Press .

Handoperated toggle action bench press ideally suitable for blanking, piercing, riveting, bending or punching, T-slotted table is provided with tool bolster for circular type die and setting plug for alignment. The ram is adjustable for tool height.

Rated capacity

Centre to back

Max. movement of handle

Clear & forward position

Clear to guide

Ram adjustment

Stroke

55 x 35 Base plate
46 prox. net weight

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003