

Product Name :
Tool Cabinet

Product Code :
NLAB-LABORATORYEM80001



Description :

Tool Cabinet

Technical Specification :

Application: Workshops Store Tools

Material: SPCC

Drawers: 1/2/3/4/5/6/7 Drawer

Usage: Garage Store Tools

Size:H760*W600*D360mmH915*W680*D460mmï¼ŒEdifferent size

Super load-bearing capacity, the whole car can load up to 500 kg

High quality ball pulley,Sliding smoothly and silently, not easy to deform, rust-proof and wear-resistant

Heavy-duty universal bearing wheel,Double bearings, silent brakes, polyurethane material, silencer, stable and wear-resistant

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003