

**Product Name :**  
Tool Kit

**Product Code :**  
NLAB-LBEMALAB230044



**Description :**

Tool Kit

**Technical Specification :**

Set of 6 Roll set of miniature open-end wrenches 3.2 to 13 mm

Set of 4 Insulated screwdrivers 3.5x75 – 4x100 -5.5x150 – 6.5x150 mm

Set of 3 Insulated Phillips screwdrivers no. 0 – 1 – 2

Set of 3 Insulated Pozidriv screwdrivers no. 0 – 1 – 2

Multi blade screwdrivers 1.5 – 2 and 2.5 mm

Tuning screwdriver for slotted heads 2 and 2.4 mm

Set of 9 Nut spinners 3.2 to 10 mm

Precision diagonal cutters

Diagonal cutters

Short-nose pliers

Angled half-round nose pliers

Crimping pliers

Soldering iron 30W

Reel of solder wire

Disordering pump

Electricians scissors

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003