

Product Name :
Tor Steel Rod Cutting Machine

Product Code :
NLAB-MACHINELABWK23017



Description :

Tor Steel Rod Cutting Machine

Technical Specification :

Heavy duty box section, M.S. Body, Double leverage type for quick cutting of light rods, lever retaining spring to hold handle failing off accidentally. and 24 are single leverage type. .

Round bars upto

mm

Square bars upto

upto

Size of sq.
blades

Approximate weight

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003