


Product Name :
Torsion Apparatus

Product Code :
NLAB-ENGINEERINGLB26002



Description :

Torsion Apparatus

Technical Specification :

Torsion Apparatus Lab Equipments Manufacturer, Suppliers & Exporters Description of Torsion Apparatus:

Torsion Apparatus is a most useful piece of bench mounted equipment and consists of a horizontal unit to take specimens up to 6 mm diameter by 450 mm long (0.25 inches x 18 inches). The torsion head is supported in ball bearings and torque is applied by a mass hanger and masses. Angular displacement is measured by a dial and pointer which can slide along the base to enable readings to be taken at any point. Supplied complete with a stainless steel specimen of 6 mm dia and one pointer.

Applications of Torsion Apparatus:

Demonstration of torsion characteristics of small diameter rods of different materials

Determination of modulus of rigidity of different materials

Features of Torsion Apparatus:

Compact bench top unit

Simple operation, with accurate repeatable results

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003