

**Product Name :**  
Traction Dynamometer

**Product Code :**  
NLAB-EDUCATIONTVET15016



**Description :**

Traction Dynamometer For Adhesion Test With Manometer of Capacity 16 Kn, Digital Display

**Technical Specification :**

Requires a drill!

Electric hammer drill

Accessories:

Metal discs for adhesion tests, diameter 20 mm, 10 pcs (consumable)

Core drill, diameter 20 mm

Adhesive glue 300 ml

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

