

Product Name :
Training System Level Control

Product Code :
NLAB-ENGINEERINGLB54007



Description :

Training System Level Control

Technical Specification :

- Experimental unit for control engineering experiments.
- Level control process with the transparent tank.
- Speed-controlled pump.
- Level measurement by the pressure sensor.
- Tank with overflow and graduated scale.
- Process schematic on the front panel.
- Disturbance variables generated by an electromagnetic proportional valve in tank outlet.

Pump

- Power consumption: 18W
- Max. flow rate: 8L/min
- Max. head: 6m
- Pressure sensor: 0...30mbar (0...300mm)
- Process schematic with controller type selection
- Proportional valve: Kvs: 0,7m³/h

Technical Data

- Level-controlled tank
- Capacity: 1200mL
- Storage tank
- Capacity: 3700mL

We are well-known manufacturers, OEM suppliers of Training System Level Control for Applied Machines. Contact us for high quality Training System Level Control for Applied Machines for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003