

**Product Name :**  
Transformer Oil Set

**Product Code :**  
NLAB-ELECTRICALAB200001



**Description :**

Transformer Oil Set

**Technical Specification :**

Transformer Oil Set Manufacturer and Supplier

This set up will consists of following facilities:

It can be operated Manual/ Motorized.

Input 1 ph 220V and output 0-60KV/100/KV with set brakedown voltag.

High voltage transformer manufactured with appoxy casted.

Oil cup provided with go and no-go gauses.

Built in control circuit and production circuit duly wired.

Built in Power on Indicator trip status indicator.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003