

Product Name :
Transformer Oil Tester

Product Code :
NLAB-TESTINGLABM23002



Description :

Transformer Oil Tester

Technical Specification :

Features :

Breakdown test of user defined

Proof test and Calibration interface

Technical Specifications :

Output voltage : AC 0 ~ 60KV / AC 0 ~ 80KV / AC 0 ~ 100KV

Rate of voltage rise : 2KV/s, $\pm 10\%$

Measuring Precision : 3 Class

Times : 1 ~ 9

Mixing Time : 0 ~ 99 second

Placing Time : 0 ~ 9 minute 59 second

Save 100 groups of data most

Voltage limit setting for proof test : 20KV ~ 70KV

Hold time for proof test : 1 ~ 10 minutes

Power source : 220 ~ 230V AC, 50Hz, 1 Phase

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003