

**Product Name :**  
Transilluminator

**Product Code :**  
NLAB-EDUSCHOOLEQ46020



**Description :**

Transilluminator

**Technical Specification :**

Being experience in our field, we are able to provide high quality Transluminators widely used in molecular biology labs.

The main application of the Transluminator is in viewing of DNA and RNA to that has been separated by electrophoresis through an agarose gel.

Premium quality equipments are used in the installation of these Transluminators.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003