

Product Name :
Transparent Water Bath

Product Code :
NLAB-SCIENTIFICEQ34003



Description :

Transparent Water Bath

Technical Specification :

Transparent Water Bath

Used for precision constant temperature and auxiliary.

Characteristics:

Using innovative design structure, effectively prevent steam leakage.

Intelligent temperature control instrument with high brightness digital display, light-touch type key has the advantages of simple operation and high precision temperature control.

Shell, interior and cover of the product are all made of transparent engineering plastics, can clearly see the object to be heated.

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com