

Product Name :
Transport Incubator

Product Code :
NLABMEDICAL-M2740007

Description :

Transport Incubator

Technical Specification :

Be it inter or intra hospital transfer, proper transportation is crucial from the time infant is picked up till reaching the destination.

Microprocessor controlled system

LED display

Double wall made of high optical quality acrylic

Precise distribution of heat over the whole hood

Removable parts easy to clean

Internal Noise Level <60 dB

Memory foam mattress for newborn comfort

Four Oval Ports

Easy access panel in front and back

Diffuser for Pipe Passage, 2 front & 2 back

Patient skin temperature probe

Passive Humidity System

Skin mode and Air control mode

Air Temperature monitoring

Patient temperature monitoring

Oxygen concentration
Comprehensive alarm system
Heating failure alarm
Skin probe failure
Air probe failure
Air circulation failure
Power supply failure
System failure
Low battery
Air micro-filter
Castors with brakes
IV Pole

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003